

MS TECH
SCIENTIFIC



Selling to
Solutioning

MS TECH SCIENTIFIC

SELLING TO SOLUTIONING

WAYEAL



MS TECH SCIENTIFIC

+91-8825969419, 7010965800

sales@mstechscientific.com

mstechscientific@gmail.com

**Primex Infra Structure Lavender Street
Kolapakkam Chennai - 600 116**


HPLC 3200

MS TECH SCIENTIFIC
WAYEAL



LC3200 series High Performance Liquid Chromatography

Main Stream split design, complete system testing function, to meet most testing needs. Customized acetonitrile-resistant check valve. The new cam curve design adopts electronic pulsation suppression technology and solvent compression to further reduce pressure fluctuation and ensure higher repeatability of test results. Standard configuration of plunger rod online cleaning function. 1800 L/mm grating imported from France. Fully Intelligent counter-control chromatography workstation



APPLICATIONS




UHPLC 3600

LC3600 series High Performance Liquid Chromatography

Innovative Ultra High Performance Liquid Chromatography Ultra-high pressure linear motor pump with 125MPa pressure resistance. Needle-in-loop sampler for loss less injection and minimal cross-contamination Liquid-core waveguide flow cell technology for lossless optical energy transmission and ultra-high detection rate Fast detection speed, high sensitivity and high separation. Flexible configuration of DAD, FLD, ELSD, RID and other detectors to meet the needs of different customers PCT international patent

WORKSTATION: SMART LAB CDS

- 21CFR Part 11 compliance
- Simple and effective data and report management
- User Friendly
- Effectively guaranteeing the safety of user data.
- High reliability for user experiments.



IC6200

MS TECH SCIENTIFIC
WAYEAL

MS TECH
SCIENTIFIC
Selling to
Solutioning

IC6200 series Ion Chromatography

The ion chromatography adopts the world's advanced technology, with reliable performance, accurate results, easy to learn and use, and intelligent. It is suitable for the routine analysis of trace and trace anions, cations, matrix complex samples and the research and development of related methods.


The IC6200 series ion chromatography, integrated design, flexible configuration, compact and convenient, the system performs high-efficiency chromatographic analysis through comprehensive application solutions and intelligent health maintenance functions, and flexible detector configurations can meet your current and future analysis needs.



IC6300

IC6300 series Intelligent Ion Chromatography

The IC6300 series intelligent ion chromatography instrument is an integrated intelligent ion chromatography system, which fully opens a new era of ion chromatography. The system is rich in configuration, powerful, simple in layout and easy to operate. The new intelligent control system, equipped with an online eluent generator, and a built-in quaternary gradient pump made of all PEEK materials, is an ideal experimental partner for you to easily achieve your experimental goals.



IC6600

MS TECH SCIENTIFIC
WAYEAL



IC6600 series Multimode Ion Chromatography

IC6600 series multi-function ion chromatography is a high-end ion chromatography system launched by Wanyi Technology in combination with comprehensive application scenarios, with a new modular design, great flexibility, more comprehensive functions, and easier operation.

By configuring conductivity detectors, amperometric detectors, and ultraviolet detectors, it can realize the detection of all items related to ion chromatography, such as conventional anions, cations, cyanide, iodine ions, sugars, small molecule organic acids, hexavalent chromium (chromic acid mist), and transition metals. The pluggable design of the amperage and conductivity detectors allows for free switching. One machine with multiple functions, while meeting the customer's routine detection, it can upgrade the functions of post-column derivation, online enrichment, online matrix elimination, etc., and its perfect and excellent performance brings chromatographic analysis to a new and higher realm. Highly flexible system to address potential challenges and advanced application scenarios, increasing productivity, expanding capacity, and improving chromatographic performance.



APPLICATIONS

AAS

AA2300 series Atomic Absorption Spectrophotometer

AA2300 series products are the fifth generation of atomic absorption spectrometer. The product is a fully automatic and intelligent spectrometer, with the continuous improvement of technology, the product automation performance and stability of the measurement to a new level, its excellent flexible configuration and rich scalability, and equipped with a powerful software system, safer and more reliable, to meet the detection needs of all walks of life.



To Know More About
www.mstechscientific.com